

AIR-FLOW® POWDERS

ORIGINAL AIR-FLOW® THERAPY
→ SUPRAGINGIVAL

CLINICAL APPLICATIONS

**AIR-FLOW®
HANDPIECES**



AIR-FLOW® > Sodium bicarbonate
POWDER > Average grain size ~65µm
CLASSIC > Available in 6 different tastes
 300 g bottle

- REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS
- SURFACE PREPARATION BEFORE DENTAL TREATMENTS



AIR-FLOW® > Sodium bicarbonate
POWDER > Average grain size ~40µm
CLASSIC
COMFORT 250 g bottle

- REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS WITH **OPTIMIZED PATIENT COMFORT**
- SURFACE PREPARATION BEFORE DENTAL TREATMENTS



AIR-FLOW® > Glycine-based
POWDER > Average grain size ~65µm
SOFT 200 g bottle

- REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS
- TOOTH POLISHING



AIR-FLOW®

→ SUB- AND SUPRAGINGIVAL



AIR-FLOW® > Erythritol-based
POWDER > Average grain size ~14µm
PLUS > Contains chlorhexidine (0,3%)
 100 g bottle

- REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS WITH **OPTIMIZED PATIENT COMFORT**
- TOOTH POLISHING
- REMOVAL OF BIOFILM FROM POCKETS **DOWN TO 4MM** ON NATURAL TEETH AND IMPLANTS
- REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, **FROM 4 TO 10MM** ON NATURAL TEETH AND IMPLANTS



**AIR-FLOW®
AIR-FLOW® PLUS**

→ SUBGINGIVAL



AIR-FLOW® > Glycine-based
POWDER > Average grain size ~ 25µm
PERIO 120 g bottle

- REMOVAL OF BIOFILM FROM POCKETS **DOWN TO 4MM** ON NATURAL TEETH AND IMPLANTS
- REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, **FROM 4 TO 10MM** ON NATURAL TEETH AND IMPLANTS



**PERIO-FLOW®
AIR-FLOW® PERIO**

NEW



NEW



"I FEEL GOOD"

EMS-SWISSQUALITY.COM

EMS
ELECTRO MEDICAL SYSTEMS SA
Chemin de la Vuarpillière 31
CH-1260 Nyon, Switzerland

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com

© EMS SA FA-510/EN Edition 08/2013

**WHAT
FOR
WHAT**

AIR-FLOW® POWDERS

EMS